

ZetaSep FPLC™ Desalting Columns

For desalting, removal of small molecules, and buffer exchange using liquid chromatography systems



ZetaSep FPLC Desalting Columns are designed for:

- Separating larger biomolecules (i.e. proteins such as antibodies, enzymes or larger nucleic acids) from unwanted smaller molecules.
- Buffer exchange, desalting, removal of low molecular weight contaminants, and reaction terminations.
- Simple, rapid and reproducible separation using a syringe, pump or liquid chromatography system.


The fractionation range for globular proteins is between 1 and 5 kD.

The molecular weight cut-off (MWCO) is approximately 5 kD, which ensures efficient separation of proteins/peptides/biomolecules larger than 5 kD from lower molecular weight molecules of less than 1 kD.

ZetaSep Desalting FPLC Columns contain **Zetadex-25 Superfine**, a beaded composite size-exclusion matrix developed and manufactured by **emp BIOTECH**. It exhibits high flow rates, excellent resolution and chemical stability. Buffer and pH effects on resolution are minimal.



ZetaSep FPLC™ Desalting Columns



High Performance Results:

Sample: 1 ml of 2 mg/ml BSA & 100 μ M of 5-Carboxyfluorescein in PBS pH 7.4 (0.05 % NaN_3).

Flow rate: 2 ml/min.

Eluent: PBS pH 7.4 (0.05 % NaN_3)

Detection: Abs. at 280 nm and 490 nm

Specifications

Column bed volume	5 ml
Size of eluted Proteins	> 5 kD
System compatibility	<ul style="list-style-type: none"> Automated liquid chromatography systems (MPLC, FPLC™, ÄKTA™, etc.) Peristaltic pump Syringe
Column dimensions	1.6 cm inner diameter x 2.5 cm height
Column body material	Polypropylene
Column ports	Inlet 10–32 (1/16") female Outlet 10–32 (1/16") male
Support	Zetadex-25 Superfine
Bead size	40–110 μ m (hydrated)
Maximum back pressure	3 bar (0.3 MPa)
Recommended flow rate	1 to 5 ml/min
Maximum recommended flow rate	10 ml/min
Storage temperature	ambient
Storage solution	20 % (v/v) ethanol

emp Biotech acknowledges and recognizes ownership of trademarks used in this publication by the respective owners of said trademarks.

Catalog Number	Description	Pack Size
ZS-0101-M005.0-005	ZetaSep FPLC Desalting	5 x 5 ml Columns
ZS-0101-M005.0-100	ZetaSep FPLC Desalting	100 x 5 ml Columns
ZS-0102-M001.0-005	ZetaSep FPLC Desalting	5 x 1 ml Columns
ZS-0102-M001.0-100	ZetaSep FPLC Desalting	100 x 1 ml Columns

