## aokin mycontrol detection kits

## Rapid and quantitative detection of mycotoxin content in food and feedstuffs

aokin mycontrol kits are for rapid and quantitative detection of mycotoxin content in food and feedstuffs. aokin mycontrol kits are designed to be used with the aokin spectrometer FP470 and the aokin Liquid Handling Workstation LHW03 for a fully automated and hands-free detection. Standardized sample preparation and mycotoxin extraction protocols are included.



ZON
AFLA
OTA
FUM
T2/HT2



## aokin mycontrol detection kits

aokin mycontrol DON Quantitative Assay for detection of Deoxynivalenol (DON)

Extraction 3 minutes \*1
Detection 3 minutes

Range of detection  $50-1250 \mu g/kg / 100-2500 \mu g/kg / 200-5000 \mu g/kg$  Commodity's Cereals (wheat, corn,...), pastries, animal feed

aokin mycontrol ZON Quantitative Assay for detection of Zearalenone (ZON)

Extraction 8 minutes \*1
Detection 3 minutes

Range of detection 25–250 μg/kg / 50–500 μg/kg / 100–1000 μg/kg Commodity's Cereals (wheat, corn,...), pastries, animal feed

aokin mycontrol AFLA Quantitative Assay for detection of Aflatoxin (Afla)

Extraction Food 25 minutes \*2
Extraction Feed 5 minutes \*1
Detection 3 minutes

Range of detection\*<sup>2</sup>  $0.1-3.5 \,\mu g/kg / 0.2-7.0 \,\mu g/kg / 0.4-14 \,\mu g/kg$ Range of detection\*<sup>1</sup>  $5-75 \,\mu g/kg / 10-150 \,\mu g/kg / 20-300 \,\mu g/kg$ Commodity's Cereals, pastries, coffee, nuts, spices, animal feed

aokin mycontrol OTA Quantitative Assay for detection of Ocratoxin (OTA)

Extraction 25 minutes \*2
Detection 3 minutes

Range of detection  $0.8-8.0 \,\mu g/kg / 1.6-16 \,\mu g/kg / 3.2-32 \,\mu g/kg$  Commodity's Cereals, pastries, coffee, nuts, spices, animal feed

aokin mycontrol FUM Quantitative Assay for detection of Fumonisin (Fum)

Extraction 15 minutes \*2
Detection 3 minutes
Range of detection 10-1250 µg/kg

Commodity's Cereals (wheat, corn,...), pastries, animal feed

aokin mycontrol T2/HT2 Quantitative Assay for detection of T2/HT2 Toxin (T2/HT2)

Extraction 8 minutes \*1
Detection 3 minutes

Range of detection  $20-400 \mu g/kg / 40-800 \mu g/kg$ 

Commodity's Cereals (wheat, corn,...), pastries, animal feed

\*¹using aokin QuickClean (SPE) columns \*²using aokin ImmunoCleanM (IAC) columns

**aokin** AG · Robert-Rössle-Str. 10 · 13125 Berlin · Germany Tel. +49 30 94 89 21 60 · Fax +49 30 94 89 21 61 · info@aokin.com

aokin is ISO 9001 certified and works according to GLP and ISO 17025.
Errors and changes reserved. aokin AG 02/2014